Optiflow™+ Nasal Cannula

- REF MYOPT9SMALL
- REF MYOPT9MEDIUM
- REF MYOPT9LARGE


Thank you for choosing F&P myAirvo™

Please read this entire guide before you first use your nasal cannula kit. Following the steps in this guide and referring to the images and diagrams will assist you with daily use of your F&P myAirvo device and its consumables.


LANGUAGES

PAGE	LANGUAGE		
1	English (en)	40	한국어 (ko)
4	Svenska (sv)	43	Norsk (no)
7	العربية (ar)	46	Nederlands (nl)
10	Čeština (cs)	49	Polski (pl)
13	Dansk (da)	52	Português do Brasil (ptbr)
16	Deutsch (de)	55	Português (pt)
19	Ελληνικά (el)	58	Română (ro)
22	Espanol (es)	61	Русский (ru)
25	Suomi (fi)	64	Slovenčina (sk)
28	Français (fr)	67	ภาษาไทย (th)
31	עברית (he)	70	Türkçe (tr)
34	Italiano (it)	73	繁體中文 (zht)
37	日本語 (ja)		


i

iii

iv

 Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

 **Australia (AU)** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号电话: +86 20 32053486 **Colombia (CO)** Tel: +57 3142852934 **France (FR)**  Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Indonesia (ID)** Tel: +62 21 509 13476 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Sweden (SE)** Fisher & Paykel Healthcare AB, Svetsarvägen 15, 2tr, 17141 Solna, Sweden Tel: +46 8 564 76 680 **Switzerland (CH)**  Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)**  Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

 **Austria (AT)** Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 45 25 70 85 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 910 38 81 18

myAirvo Optiflow+ Nasal Cannula

[REF] MYOPT9SMALL
Optiflow+ Nasal Cannula Small

[REF] MYOPT9MEDIUM
Optiflow+ Nasal Cannula Medium

[REF] MYOPT9LARGE
Optiflow+ Nasal Cannula Large

INDICATIONS FOR USE

The F&P Optiflow+ is a nasal cannula patient interface for use with specified respiratory gas humidifiers to treat spontaneously breathing pediatric (3 years and older) to adult patients who would benefit from receiving high-flow warmed and humidified respiratory gases to the upper airway.

This device is designed to be used in a long-term care facility by appropriately qualified healthcare professionals, or in the home by lay users operating the device as prescribed by a healthcare professional.

SPECIFICATIONS

- myAirvo Series humidifiers with MYAIRVOKIT1/MYAIRSPIRAL/900PT560/900PT560E series tubes or kits

Note: Not all products are available in all countries.

- Flow range: MYOPT9SMALL 10–60* L/min
MYOPT9MEDIUM 10–60* L/min
MYOPT9LARGE 10–60* L/min

* Maximum flow rate depends on cannula size and humidifier model. Do not adjust flow rates without first consulting your healthcare provider. Refer to the myAirvo Series User Manual for more information.

Note: Flow rates are expressed in body temperature, pressure, saturated (BTPS).

Additional information:

- This product is intended to be used for a maximum of 30 days provided that daily and weekly care instructions are followed.
- Replace the interface every 30 days, or earlier if damage or deterioration occurs.
- Ambient operating temperature: 18 °C to 28 °C (64 °F to 82 °F)
- Permissible transport and storage conditions: -10 °C to 50 °C (14 °F to 122 °F)
- Refer to the myAirvo Series User Manual for permissible operating conditions.

i SET-UP INSTRUCTIONS

Step 1

Select the appropriate size nasal cannula. The prongs

must not create a seal in the nares. A clear gap must be visible around each prong.

Step 2

Adjust the headstrap to fit. Do not over-tighten.

Step 3

- Ensure the headstrap clip is attached to the headstrap, to prevent cannula from being pulled out of the nares.
- The cannula can become unattached if not used with the headstrap clip.

Step 4

Attach the breathing tube clip to a secure location (e.g. clothing or bedding) to prevent the cannula from pulling on the ears or off the face.

ii DURING USE

- As appropriate for the patient's condition, monitor for any disruptions to the patient receiving flow.
- To prevent condensate, remove/minimize the impact of anything that may cool the tubing (e.g., a fan, air conditioning).
- Notice to user:** If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

iii DAILY CARE

Step 1

Disconnect the patient interface from the heated breathing tube.

Step 2

Rinse the interface in drinking-quality water then reconnect to the breathing tube whilst still in Drying Mode (refer to the myAirvo Series User Manual) to dry the interface.

iv WEEKLY CARE

Step 1

Handwash the interface in a tub of warm water with a mild dishwashing detergent. Ensure the interface is visually clean.

Step 2

Rinse the interface thoroughly in drinking-quality water. Ensure that all soap residue has been removed. Leave the interface to dry away from direct sunlight or heat, then reconnect it to the breathing tube.

SIDE EFFECTS

- Skin reddening may occur during nasal high flow therapy.

DISPOSAL

Dispose of products according to local guidelines. Packaging and User Manual should follow recycling guidelines as applicable. A used product should be disposed of as contaminated waste. In the home, place the product in the general waste collection. In long-term care facilities, follow the applicable waste disposal procedure.

WARNINGS

- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption to therapy, to avoid serious injury or death.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Failure to use the set-up described above can compromise performance and affect patient safety.
- If you are experiencing severe symptoms, contact your healthcare provider directly.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

**Exposure to oxygen increases the risk of fire.
To avoid risks of burns, patient injury or death:**

- Do not use system near a naked flame or any ignition source, including while smoking.
- Only use lotions and/or salves that are labeled as oxygen compatible.
- Do not lubricate connections or tubing.
- Do not leave the interface on bed coverings or cushions if the device is turned on but not in use. Oxygen can make the materials more flammable.

CAUTIONS

- For single patient use only. Reuse between patients may result in transmission of infectious substances.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.

- Do not sterilize or disinfect this product, or exceed the maximum duration of use to avoid degradation of materials and patient harm.
- Do not use if individual interface packaging is not sealed to avoid contamination.
- Deviations from the recommended methods of daily and weekly care may result in inadequate cleaning which may lead to adverse health effects.
- Using the cannula may cause skin reddening. Monitor for signs of skin damage or irritation and adjust the cannula if necessary.

Failure to comply with the following may result in loss of therapy:

- Do not crush or stretch tube. Do not obstruct or modify any part of the device.
- Do not use if tubing is damaged with holes, tears or kinks.
- Check for condensate regularly. Drain condensate away from patient as required, to prevent blockages or condensate from reaching the patient.
- Using the device for longer than the specified time can result in degradation of materials or infection.

SYMBOL DEFINITIONS

	Catalog number
	Manufacturer
	Importer
	Distributor
	Authorized European Representative
	Authorized Swiss Representative
	Authorized UK Representative
	Batch code
	Date and country of manufacture CC NZ: New Zealand MX: Mexico
	Use-by date
	Storage and transport temperature range

Rx Only	Prescription only
	Single patient multiple use
	Polyethylene Terephthalate
	Polyethylene Terephthalate, Low Density Polyethylene
	Refer to instructions
	Correct FITTING INSTRUCTION method
	Incorrect FITTING INSTRUCTION method
	Small nasal cannula size
	Medium nasal cannula size
	Large nasal cannula size
	Medical Device
CE 0123	European Conformity - TÜV SÜD
	Paper, recyclable

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