

Product data sheet

Alvera ThermaAir 10, Heated CPAP Tubing for Air 10 Series Machines | REF 55110 | GTIN 5065027735385

Technical and commercial reference for Alvera ThermaAir 10. Fields marked to be confirmed are pending final production data and will be issued in the next revision.

Identification

Brand	Alvera
Product	ThermaAir 10 heated CPAP tubing
Reference	55110
GTIN	5065027735385
Product type	Compatible heated tube, third party accessory
Placed on the market by	Elarix Medical Ltd, Unit 3, Twenty4 Secure Storage, Leeds LS12 4HP, United Kingdom
Country of manufacture	To be confirmed

Description

A heated tube for Air 10 series therapy machines. A heating element runs the length of the hose and is controlled by the machine's own climate control function, holding the delivered air above its dew point so humidity reaches the airway rather than condensing inside the tube. Supplied as a single non-sterile tube for single wearer use.

Technical data

Machine end	Integrated electrical and air connector
Mask end	Standard 22 mm cuff
Length	To be confirmed
Internal diameter	To be confirmed
Hose material	To be confirmed
Element supply	Drawn from the machine through the connector
Temperature control	Set on the machine, through its climate control function
Latex	To be confirmed
Packed weight	To be confirmed

Compatibility

Fits	Air 10 series machines, including AirSense 10, AirCurve 10 and Lumis series
-------------	---

Requires	A heated humidifier and a climate control function on the machine
Does not fit	AirSense 11, AirCurve 11 and AirMini machines

Use and life

Intended use	Single wearer, not reissued between people
Cleaning	Weekly wash in warm water with a mild detergent, connector kept dry
Renewal interval	Typically every three to six months
Storage	Loosely coiled, dry, away from direct sunlight and heat

Commercial

Pack	Single tube
Case quantity	To be confirmed
Case dimensions and weight	To be confirmed
Case GTIN	To be confirmed once the case configuration is fixed

Issued by Elarix Medical Ltd. This data sheet is provided for information. The instructions supplied with the therapy machine and with the tube take precedence.

Elarix Medical Ltd, Unit 3, Twenty4 Secure Storage, Leeds LS12 4HP, United Kingdom. Company No. 16602969.
support@elarixmedical.co.uk. www.elarixmedical.com