

prisma VENT40

Invasive and non-invasive ventilation.

With prisma VENT40 you'll have mobility and convenience at every level:

- pressure range up to 40 hPa
- use of leakage circuit system (15 mm / 22 mm)
- with autoS/T mode (autoEPAP, autoF) and target volume (100 – 2,000 ml)
- includes helpful functions for COPD treatment: AirTrap Control, maximum setting for expiratory ramp and Trigger Lockout
- integrated port for oxygen feed
- three ventilation programs stored in device memory
- complete alarm management
- rechargeable battery with 10-hour capacity
- easy to use thanks to intuitive menus, simple operation and quick access

Quiet as a whisper – for the comfort and well-being of patients and partners.

Our Software Solutions



- **prisma CLOUD**
The telemedical connection for experts and patients. Information updated daily on therapy, remote setting of therapy parameters. Fully automated, from anywhere.
- **prismaTS**
Therapy software for experts. For detailed data analysis and remote settings.





prisma VENT40

Accessories

Humidifier			
prismaAQUA	WM 29490	prismaHYBERNITE	
		breathing tube, heated	
Leakage circuit		22 mm Ø	WM 29067
22 mm Ø, black	WM 23962	15 mm Ø	WM 29083
22 mm Ø, autoclavable	WM 24667		
22 mm Ø	WM 24445	Inverter	
15 mm Ø	WM 29988	12 V	WM 24616

You'll find the complete line of accessories in the prisma VENT accessories catalog.

Technical specifications

Product class as per MDR (EU) 2017/745	IIa	Ti/Ti max Ti min, Ti max, Ti timed	0.5 s – 4 s 0.2 s – 4 s auto (only Ti timed)
Dimensions in cm (W x H x D):	21.8 x 17.5 x 21.8	Respiratory rate	5 – 60 breaths/min
Weight	2.4 kg	Target volume	100 ml – 2000 ml
Weight of internal rechargeable battery (for variant with battery)	0.63 kg	Trigger levels - Inspiration - Expiration	1 (high sensitivity) to 8 (low sensitivity) 95% to 5% of maximum flow in steps of 5%
Maximum air flow at 20 hPa	> 220 l/min	Pressure increase speed	Level 1: 100 hPa/s Level 2: 80 hPa/s Level 3: 50 hPa/s Level 4: 20 hPa/s
Operating time of internal rechargeable battery with the following settings: T mode, f = 20 /min, Ti = 1 s, PEEP = 4 hPa, Vt = 800 ml Passive lungs: Resistance R = 5 hPa (l/s); Compliance C = 50 ml/hPa	> 10 hours	Pressure decrease speed	Level 1: 100 hPa/s Level 2: 80 hPa/s Level 3: 50 hPa/s Level 4: 20 hPa/s max: slow pressure adjustment
IPAP pressure range	4 – 40 hPa ±1.2 hPa (±8% of setting)	Maximum allowed flow with oxygen supply	15 l/min
PEEP pressure range	4 – 25 hPa ±1.2 hPa (±8% of setting)		

Device designation	Article no.
prisma VENT40	WM 29500-1110
prisma VENT40 with integrated rechargeable battery	WM 29550-1110

Please refer to the currently valid price list.

CE 0197

© Protected by copyright. Reproduction of any kind requires the express permission of Löwenstein Medical. All information is subject to change for technical reasons.



Made in
Germany

Löwenstein Medical Technology
Kronsaalsweg 40
22525 Hamburg, Germany
T. +49 40 54702-0
F. +49 40 54702-461
info@loewensteinmedical.com
loewensteinmedical.com

Sales + Service
Löwenstein Medical
Arzbacher Straße 80
56130 Bad Ems, Germany
T. +49 2603 9600-0
F. +49 2603 9600-50
info@loewensteinmedical.com
loewensteinmedical.com



With people in mind